

Eyewash Station

Product Images

Product Code: S26-7238



Short Description

Comprehensive eye care first aid kit containing sterile eyewash solution which can be used to irrigate the eye or cleanse a wound.

Each 500ml bottle contains the active ingredient of 0.9% w/v sodium chloride EP and is simple to use by twisting the cap and removing the seal.

This particular kit also has a handy inbuilt mirror assisting when applying the sterile eyewash solution and also comes complete with two boxed eye dressings which can be used to protect the eye after irrigation.

The eyewash kit shows the contents at a glance and can be wall mounted for high visibility and easy access, for use in the workplace.

Description

Comprehensive eye care first aid kit containing sterile eyewash solution which can be used to irrigate the eye or cleanse a wound.

Each 500ml bottle contains the active ingredient of 0.9% w/v sodium chloride EP and is simple to use by twisting the cap and removing the seal.

This particular kit also has a handy inbuilt mirror assisting when applying the sterile eyewash solution and also comes complete with two boxed eye dressings which can be used to protect the eye after irrigation.

The eyewash kit shows the contents at a glance and can be wall mounted for high visibility and easy access, for use in the workplace.